Supercritical Fluid Application Notes

SFE543: Enzyme-Assisted Supercritical Carbon Dioxide Extraction of Piperine in Black Pepper Oleoresin

Black pepper (Piper nigrum L.) is known as the King of Spices and is the most popular spice internationally. The active ingredient of black pepper is piperine, a potent therapeutic compound.

Piperine is typically extracted by methods that may promote thermal degradation, hydrolysis of useful pharmaceutical compounds, and traces of solvents in the extracts. Supercritical CO₂ extraction offers an alternative, green extraction technique suited for solvent free extracts of spices for therapeutic applications.

In this work, enzyme-assisted supercritical carbon dioxide (SC-CO₂) extraction of black pepper oleoresin was investigated using a-amylase (from Bacillus licheniformis) for enhanced yield of piperine-rich extract possessing a good combination of phytochemical properties.

Black pepper contains approximately 58% carbohydrate and starch (30.4 %) is the predominant carbohydrate in the sample. Since starch may inhibit the extraction of piperine, the use of starch degrading enzymes, such as a-amylase for pre-treatment of the pepper matrix prior to extraction renders extraction easier and improves the piperine yield of the extracts.



Materials

(Please refer to references for detailed procedures)

Malabar Garbled black pepper

Standard piperine (97% pure)

a-amylase from B. licheniformis (lyophilized powder)

Sample Preparation

Black pepper berries were ground using an electric mixer to dp of 0.42

For enzyme-assisted extraction of black pepper oleoresin, the powdered pepper sample was mixed with the lyophilized enzyme in optimized ratio (enzyme: black pepper powder 1: 5000) and subjected to continuous modes of extraction.

Equipment

Applied Separations Supercritical Extraction Equipment: *Spe-ed* SFE 2



930 Hamilton Street Allentown, PA 18101 610-770-0900 www.appliedseparations.com

Supercritical Fluid Application Notes

Extraction Conditions

Extraction vessel: 50 mL Sample: 20 g Pressure: 300 Bar Temperature: 60 ºC CO₂ Flow Rate: 2 L/min (gas) Static time: 30 minutes Dynamic time: 15 minutes Collection: 60 mL pre-weighed vial

Analysis:

Densitometric analyses of black pepper extracts were conducted to estimate the total piperine content using a Camag TLC. The plates were scanned with Camag TLC Scanner 3 at 337 nm using a Deuterium lamp at a scanning speed of 20 mm/s. Amount of piperine present in the extracts was determined from the standard curve prepared using pure piperine.

Results

Enzyme-assisted SC-CO₂ vs SC-CO₂ extraction of black pepper at 60C and 300 bar.

Method	Yield piperine
	(g extract/100 g dry
	black pepper)
Enzyme Assisted SC-	
CO ₂	5.3 (SD .2)
SC-CO ₂	4.6 (SD .2)

Conclusion

Alpha-amylase assisted SC-CO₂ extraction of black pepper oleoresin enhanced the yield of piperine-rich extract with a good mix of phytochemical properties (antioxidant activity, total phenolic content, reducing power, and anti-inflammatory activity).

References

Dutta S, Bhattacharjee P. Enzyme-assisted supercritical carbon dioxide extraction of black pepper oleoresin for enhanced yield of piperinerich extract. Journal of bioscience and bioengineering. 2015 Jul 1;120(1):17-23.

https://pubmed.ncbi.nlm.nih.gov/25617183/

Dutta S, Bhattacharjee P. Nanoliposomal encapsulates of piperine-rich black pepper extract obtained by enzyme-assisted supercritical carbon dioxide extraction. Journal of Food Engineering. 2017 May 1;201:49-56.

Nagavekar N, Singhal RS. Enhanced extraction of oleoresin from Piper nigrum by supercritical carbon dioxide using ethanol as a co-solvent and its bioactivity profile. Journal of Food Process Engineering. 2018 Feb;41(1):e12670.

Dutta S, Bhattacharjee P, Bhattacharyya N. Assessment of shelf lives of black pepper and small cardamom cookies by metal oxide-based electronic nose using spoilage index. Food and Bioprocess Technology. 2017 Nov;10:2023-33.

Dutta S. Enzyme-assisted supercritical carbon dioxide extraction of active principles of black pepper and small cardamom and applications of the extracts for design of nutraceutical foods and supplements.



930 Hamilton Street Allentown, PA 18101 610-770-0900 www.appliedseparations.com